

Partner for
disinfection

Rw




Linotox[®]

Linen Disinfectant Concentrate

(Chlorine-free)

Linotox[®] is designed to give a complete and efficacious, Linen Disinfection Cycle. It cleanses, disinfects, all hospital linen and leaves behind a pleasant smell, removing all odours which are almost always associated with hospital linen.



Bacteria

Fungi



How
seriously
do you treat
Infected Linen?

Linotox®

Linen Disinfectant Concentrate

(Chlorine-free)



Linotox® is the perfect “pre-soak” product for contaminated linen, as well as a “final rinse” product which cleanses, and disinfects hospital linen. Its special formula contains a **leave-behind fragrance** which is pleasant, and eliminates bad odours.

Climatic conditions of the laundry section are moist and warm which contributes to the recontamination of already treated laundry, therefore, sometimes it is important to use a good disinfectant in the “final rinse”



Linotox® is ideal for :

- Infected and contaminated linen.
- Foul smelling linen with excretion and urine.
- Linen soiled with blood and sputum.
- Clothing of patients of communicable diseases.
- Reusable Surgical Dressings.
- Blankets, gowns, undergarments & floor mops.

Linotox® Offers

- Extra high cleansing power.
- Fabric protectors.
- Excellent leave-behind fragrance, and removes all bad odours.
- Chlorine-free formula, therefore, excellent for white and coloured linen.
- Blood dissolving capacity.
- Total disinfection without addition of any other cleansing agent.
- Removal of fresh stains.
- Excellent moistening capacity.

Directions for Use :

- Infected Linen** (soiled with blood, mucous, urine, stool, etc.) **squeeze out contaminants**, then PRE-SOAK in **Linotox®** 1% solution for half hour. (i.e. 10 caps full i.e. 100 ml, in 10 Litres of water).
- Non Infected Linen** (Blankets, gowns, covers, masks, etc.) PRE-SOAK in 0.5% **Linotox®** solution, for half hour (5 caps full i.e. 50 ml. in 10 Litres of water).

Linotox® ensures efficacy against :

- All Bacteria
- All Fungi

NOTE :

- Please do not use any other cleansing agent.
- Contains fabric protectors, therefore, safe on all linen.
- Can be used in washing machine cycle also.
- Can be used in the “Final-rinse” for best results.
- Safe for coloured linen also.

Composition :

Each 100 gms contain :

Didecyl Dimethyl Ammonium Chloride	5.0 gms
Perfume	0.5 gms

Packing :

1 Bottle. :	500 ml.
1 Carton :	20 bottles.

